


BREAST PUMP SET-UP INSTRUCTIONS

1. A breast pump should be in your room. If there is not one, please check with your nurse.
2. The membranes and caps are placed over the pumping mechanism.



3. Connect the tubing to the breast pump. One end goes to the caps on top of the pump and the other goes to the flanges of the milk collection set-up.



4. Make certain that the white membrane disk is intact and seated in the yellow valve. The yellow valve connects to the flanges of the milk collection set-up.



5. Attach the milk collection containers to the flanges. Turn the pump on, adjust the suction to a level of comfort and pump as directed by the lactation nurse.

